





PHLEBOLOGY

LEONARDO® DEPTH CONTROL



Technology
meets Anatomy



High combined power and energy
in the junction area
for optimal sealing


Adjustment of **both wavelengths-parameters**
for extrafascial side branches

Adaptation to popliteal fossa region with regard to
muscle veins and nerves


Perforator treatment
with adjusted
higher energies

Identification of nerves and adaptation
to the small saphenous vein
target parameters

The art of water absorption



Combined Efficacy – 1940 nm and 1470 nm joining together



Simplified action mechanism

Araujo et al Comparison of the effects of endovenous laser ablation at 1470 nm versus 1940 nm and different energy densities DOI:10.1177/0268355518778488

1470 nm – Proven and effective

- Gold standard for safe & durable vein closure
- Also for large vein diameters (>15 mm)

1940 nm – Greater precision, gentler touch

- Lower penetration depth than 1470 nm – optimal control
- Particularly suitable for sensitive anatomical areas


Blending – It's the mix

Combining both wavelengths enables precise tailoring of the laser–tissue interaction to optimize recovery time, minimize side effects, and ensure reliable vein closure.

The energy delivery can be adapted to each vein's diameter, wall thickness, and proximity to nerves or other sensitive structures, allowing controlled thermal confinement. This dual-wavelength approach effectively limits lateral heat diffusion into surrounding tissue, thereby further reducing pain, hematomas, and nerve irritation.

Safety and comfort for your patients

- Less post-operative discomfort
- Faster recovery and less downtime
- Ultimate control





LEONARDO®

Model	LEONARDO® DUAL 20
REF	LL1470+1940nm20W
Wavelength	1470 nm ± 30 nm + 1940 nm ± 30 nm
Power	20 Watt (12 W 1470 nm + 8 W 1940 nm)
Fiber diameter	≥ 360 µm
Target beam	650 nm, red 4 mW, user controlled intensity
Treatment mode	CW, Pulse Mode, ELVeS® Signal, ELVeS® Segment
Pulse length / - break	0.01 – 60 sec / 0.01 – 60 sec
Power supply	100 – 240 VAC, 50/60 Hz, 450 VA
Dimensions (H x W x D)	approx. 28 cm x 37 cm x 9 cm
Weight	approx. 8.5 kg

All laser sets include three safety goggles, a foot switch, a door interlock, a power cable, and a manual in a carrying case.



Technical highlights

- Power: Up to 12 W (1470 nm) / 8 W (1940 nm)
- Modes: Single wavelength or combined use
- Display: Intuitive user interface with touchscreen
- Portable: Compact, lightweight design – ideal for clinics and practices

biolitec® worldwide

biolitec biomedical technology GmbH

Jena, Germany
Phone: +49 3641 519 53 0

biolitec Schweiz GmbH

Wollerau, Switzerland
Phone: +41 55 555 30 20

biolitec España SL

Madrid, Spain
Phone: +34 91 9910857

biolitec Italia SRL

Milano, Italy
Phone: +39 02 8423 0633

biolitec Tıbbi Cihazları Ltd. Şti.

Istanbul, Turkey
Phone: +90 216 574 7456

OOO biolitec Spb

Saint-Petersburg, Russia
Phone: +7 812 4493752

biolitec FZ LLC

Dubai, UAE
Phone: +971 44 29 85 92

biolitec laser science and technology Shanghai Ltd.

Shanghai, China
Phone: +86 21 6308 8856

biolitec Sdn. Bhd.

Selangor, Malaysia
Phone: +60 3 5569 7158

biolitec India Private Ltd.

Bangalore, India
Phone: +91 93275 11005

PT. Biolitec

Tangerang, Indonesia
Phone: +62 21 295 57 419

biolitec Korea Ltd.

Seoul, Republic of Korea
Phone: +82 2 701 4707

Equipos Laser de Uso Médico y Fibra Óptica SA de CV

(Biolitec Mexico & Latín América)
Phone: +52 (55) 5573 1800

Biolitec Biotecnologia Comércio, Importação, Exportação LTDA

São Paulo, Brazil
Phone: +55 11 2093 8602



biolitec GmbH

Untere Viaduktgasse 6/9
A-1030 Wien
info@biolitec.com

Legal Manufacturer

CeramOptec GmbH
Siemensstr. 44
53121 Bonn

Find your
local biolitec®
partner

