



LEONARDO® Mini Athletic Laser



Technical specifications	LEONARDO® Mini Athletic
REF	SL810+980nm15W
Power/Wavelength	7.5 Watt (810 nm) / 7.5 Watt (980 nm)
Fiber diameter	≥ 360 µm
Aiming beam	635 nm, max. 4 mW
Treatment mode	CW, Pulse Mode (optional)
Pulse duration/-break	0.01 - 60 sec. / 0.01 - 60 sec.
Power supply	110 - 240 VAC, 50 - 60 Hz (7.2 VDC @ 36 Watt)
Batteries	Li-ion Batteries
Dimensions (H × W × D)	6 cm × 9 cm × 21.5 cm
Weight	900 g

Fiber

REF	Produkt	Qty *	Length [m]	OD ø [mm]
503201921	Reusable Bare Fiber 600 µm, Flat Tip, IC (10 × 1 h)	5	3	0.95

* Packaging unit

