

UROLOGY

LEONARDO[®]
DUSTER

Fibers and Accessories



Technology
meets Anatomy

Leonardo® Duster ScopeSafe Fiber- Single Use

REF	Product	Length [m]	Core ø [µm]	OD ø [µm]	Connector	NA
503200100	Leonardo® Duster ScopeSafe Fiber 150µm, IC	3.1	150	315	IC	0.26
503200105	Leonardo® Duster ScopeSafe Fiber 200µm, IC	3.1	200	400	IC	0.28
503200110	Leonardo® Duster ScopeSafe Fiber 272µm, IC	3.1	272	460	IC	0.22
503200115	Leonardo® Duster ScopeSafe Fiber 365µm, IC	3.1	365	550	IC	0.22
503200120	Leonardo® Duster ScopeSafe Fiber 550µm, IC	3.1	550	750	IC	0.22
503200125	Leonardo® Duster ScopeSafe Fiber 800µm, IC	3.1	800	1200	IC	0.22
503200130	Leonardo® Duster ScopeSafe Fiber 1000µm, IC	3.1	1000	1400	IC	0.22

Leonardo® Duster Reusable Fibers

REF	Product	Length [m]	Core ø [µm]	OD ø [µm]	Connector	NA
503200150	Leonardo® Duster Fiber 200µm, Reusable, IC	3.2	200	400	IC	0.28
503200155	Leonardo® Duster Fiber 272µm, Reusable, IC	3.2	272	460	IC	0.22
503200160	Leonardo® Duster Fiber 365µm, Reusable, IC	3.2	365	550	IC	0.22
503200165	Leonardo® Duster Fiber 550µm, Reusable, IC	3.2	550	750	IC	0.22
503200170	Leonardo® Duster Fiber 800µm, Reusable, IC	3.2	800	1200	IC	0.22
503200175	Leonardo® Duster Fiber 1000µm, Reusable, IC	3.2	1000	1400	IC	0.22

The IC abbreviation in the product name stands for "Inter-Company", which indicates an IC-connector that can only be used with Leonardo® laser series of the Biolitec group.
All Leonardo® Duster reusable fibers can be used up to 10 times.

Accessories for the Reusable Fibers

REF	Product	Core ø [µm]
400100104	Medi Strip 0.3/0.6 BF 150µm, autoclavable	200
400100105	Medi Strip Bare Fiber 200µm, autoclavable	272
400100110	Medi Strip Bare Fiber 400µm, autoclavable	365
400100115	Medi Strip 0.7/1.2 BF 600µm, autoclavable	550
400100116	Medi Strip 0.95/1.5 BF 800µm, autoclavable	800
400100120	Medi Strip 1.0/1.5 BF 1000µm, autoclavable	1000
400100130	Ceramic Fiber Cleaver	All Fibers



biolitec Holding GmbH & Co KG
Untere Viaduktgasse 6/9
A-1030 Wien
Phone +43 1 3619 909 50
info@biolitec.com

Legal Manufacturer
CeramOptec GmbH
Siemensstr. 44
53121 Bonn

Imprint
biolitec Holding GmbH & Co KG
Untere Viaduktgasse 6/9
A-1030 Wien
Phone: +43 1 3619 909 50
www.biolitec.com

All fibers are to be used only with the
1940 nm Leonardo® Duster Laser.