

UROLOGY

LEONARDO[®]
DUSTER

Fibers



Technology
meets Anatomy

Leonardo® Duster ScopeSafe Fiber

REF	Product	Length [m]	Core ø [μm]	OD ø [μm]	Connector	NA
503200100	Leonardo® Duster ScopeSafe Fiber 150μm, IC	3.1	150	315	IC	0.26
503200105	Leonardo® Duster ScopeSafe Fiber 200μm, IC	3.1	200	400	IC	0.28
503200110	Leonardo® Duster ScopeSafe Fiber 272μm, IC	3.1	272	460	IC	0.22
503200115	Leonardo® Duster ScopeSafe Fiber 365μm, IC	3.1	365	550	IC	0.22
503200120	Leonardo® Duster ScopeSafe Fiber 550μm, IC	3.1	550	750	IC	0.22
503200125	Leonardo® Duster ScopeSafe Fiber 800μm, IC	3.1	800	1200	IC	0.22
503200130	Leonardo® Duster ScopeSafe Fiber 1000μm, IC	3.1	1000	1400	IC	0.22

The IC abbreviation in the product name stands for "Inter-Company", which indicates an IC-connector that can only be used with Leonardo® laser series of the Biolitec group.



biolitec Holding GmbH & Co KG
Untere Viaduktgasse 6/9
A-1030 Wien
Phone +43 1 3619 909 50
info@biolitec.com

Legal Manufacturer
CeramOptec GmbH
Siemensstr. 44
53121 Bonn

Imprint
biolitec Holding GmbH & Co KG
Untere Viaduktgasse 6/9
A-1030 Wien
Phone: +43 1 3619 909 50
www.biolitec.com

All fibers are to be used only with the
1940 nm Leonardo® Duster Laser.