


# ELVeS<sup>®</sup> smart

## Pullback device



Continuous and homogenous withdrawal of the ELVeS<sup>®</sup> Radial fiber family during the endovenous laser ablation of varicose veins.

ELVeS<sup>®</sup> smart, your add-on for a perfect match using ELVeS<sup>®</sup> Radial fibers and our reliable biolitec<sup>®</sup> diode lasers with ELVeS<sup>®</sup> Signal Mode. This leads to unrivalled homogeneity of your treatment, all this with the push of a button.



## ELVeS<sup>®</sup> smart

REF	L-ELVeSmart-Set
Pulling speed	0.5; 1.0; 1.5; 2.0 mm/s
Fibers	ELVeS Radial <sup>®</sup> , ELVeS Radial 2ring <sup>®</sup> ELVeS Radial <sup>®</sup> slim, ELVeS Radial 2ring slim
Power Supply	100 – 240 VAC, 50 – 60 Hz
Dimensions (H × W × D)	12 cm × 16.5 cm × 22 cm
Weight	2.4 kg

