


ELVeS[®] smart

Pullback device



Continuous and homogenous withdrawal of the ELVeS[®] Radial fiber family during the endovenous laser ablation of varicose veins.

ELVeS[®] smart, your add-on for a perfect match using ELVeS[®] Radial fibers and our reliable biolitec[®] diode lasers with ELVeS[®] Signal Mode.



ELVeS[®] smart

| | |
|------------------------|---|
| REF | L-ELVesSmart-Set |
| Pulling speed | 0.5; 1.0; 1.5; 2.0 mm/s |
| Fibers | ELVeS Radial [®] , ELVeS Radial 2ring [®] ELVeS Radial [®] slim, ELVeS Radial [®] 2ring slim |
| Power Supply | 100 - 240 VAC, 50 - 60 Hz |
| Dimensions (H x W x D) | 12 cm x 16,5 cm x 22 cm |
| Weight | 2,4 kg |