


ELVeS[®] smart

Pullback device



Continuous and homogenous withdrawal of the ELVeS[®] Radial fiber family during the endovenous laser ablation of varicose veins.

ELVeS[®] smart, your add-on for a perfect match using ELVeS[®] Radial fibers and our reliable biolitec[®] diode lasers with ELVeS[®] Signal Mode. This leads to unrivalled homogeneity of your treatment, all this with the push of a button.



ELVeS[®] smart

REF	L-ELVeSsmart-Set
Pulling speed	0.5; 1.0; 1.5; 2.0 mm/s
Fibers	ELVeS Radial [®] , ELVeS Radial 2ring [®] ELVeS Radial [®] slim, ELVeS Radial 2ring slim
Power Supply	100 – 240 VAC, 50 – 60 Hz
Dimensions (H × W × D)	12 cm × 16.5 cm × 22 cm
Weight	2.4 kg

