


Our products



LEONARDO®



Model	LEONARDO® Mini 1470 nm	LEONARDO® Mini Dual	LEONARDO® DUAL 45
REF	SL1470nm12W	SL980 + 1470nm14W	SL980 + 1470nm45W
Wavelength	1470 nm	980 nm and 1470 nm	980 nm and 1470 nm
Power	12 W (1470 nm)	10 W (980 nm) / 4 W (1470 nm)	max. 45 Watt (1470 nm / 15 Watt + 980 nm / 30 Watt) separately adjustable
Fiber diameter	≥ 360 µm	≥ 360 µm	≥ 360 µm
Aiming beam	635 nm, max. 4 mW	635 nm, max. 4 mW	532 nm and 635 nm, green 1 mW, red 4 mW, user controlled intensity
Treatment mode	CW, Pulse Mode (optional), ELVeS® Signal	CW, Pulse Mode (optional)	CW, Pulse Mode, ELVeS® Signal, ELVeS® Segment, Derma Mode
Pulse duration /-break	0.01 – 180 sec. / 0.01 – 180 sec.	0.01 – 180 sec. / 0.01 – 180 sec.	0.01 – 60 sec / 0.01 – 60 sec
Power supply	110 – 240 VAC, 50 - 60 Hz (12 VDC @ 100 W)	110 – 240 VAC, 50 – 60 Hz (12 VDC @ 65 W)	110 – 240 VAC, 50 / 60 Hz, 450 VA
Batteries	Li-ion batteries	Li-ion batteries	–
Dimensions (H × W × D)	6.0 cm × 9.0 cm × 21.5 cm	6.0 cm × 9.0 cm × 21.5 cm	approx. 28 cm × 37 cm × 9 cm
Weight	900 g	900 g	approx. 8.5 kg

All laser sets incl. 3 safety goggles, foot switch, interlock connector, power cord and manual in a carrying case.

Fibers



REF	Product	PU [Packaging unit]	length [m]	ø fiber [mm]
503100250	FiLaC® Fistula Probe, IC	10	2.6	1.85
503200740	Bare Fiber 600 µm, Flat Tip, IC	10	2.6	0.96

Kits



503100220	LHP® Procedure Kit, IC	5	2.6	1.85
503100255	FiLaC® Fistula Kit, IC	5	2.6	1.85

Accessories

REF	Product	PU [Packaging unit]
400100100	Universal Dual Luer Handpiece	1
LA1371	Laser Safety goggles 950 – 1010 L4 + 1470 L2 (FULL), transparent	1

