


# Holmium Fibers

Outstanding high-quality  
Holmium laser fibers




- Exceptional quality
- High stability and durability
- Suitable for various laser systems



## biolitec Holmium laser fibers

With many years of expertise, biolitec manufactures these premium fibers in own production. Your clinic and your patients profit from our high-tech fibers of exceptional quality.



### Single-use

REF	Product	Length [m]	Core $\varnothing$ [ $\mu$ m]/[Fr]	OD $\varnothing$ [ $\mu$ m]/[Fr]	Connector	NA
500200326	Megabeam HBF5F 230 - 253, OD 420 $\mu$ m	3.1	272	420	SMA905	0.28
500200327	Megabeam HBF5F 365 - 403, OD 550 $\mu$ m	3.1	365	550	SMA905	0.22/0.37
500200328	Megabeam HBF5F 550 - 603, OD 750 $\mu$ m	3.1	550	750	SMA905	0.22/0.37

### Re-use

REF	Product	Length [m]	Core $\varnothing$ [ $\mu$ m]/[Fr]	OD $\varnothing$ [ $\mu$ m]/[Fr]	Connector	NA
500200789	Megabeam RHB5F 230/253 - 3	3.0	272	420	SMA905	0.22/0.36
500200790	Megabeam RHB5F 365/400 - 3	3.0	365	550	SMA905	0.22/0.36
500200791	Megabeam RHB5F 550/605 - 3	3.0	550	750	SMA905	0.22/0.36
500200792	Megabeam RHB5F 945/1000 - 3	3.0	945	1400	SMA905	0.22/0.36

### Side Fiber

REF	Product	Length [m]	Core $\varnothing$ [ $\mu$ m]/[Fr]	OD $\varnothing$ [ $\mu$ m]/[Fr]	Connector	NA
500300451	Megabeam SF - 2100 - H	3.2	550	950	SMA905	0.22

fibers are compatible with, but not limited to, the StoneLight®, COHERENT®, Deka, LUMENIS®, Dornier MedTech and NewStarLaser systems.