

Our products






LEONARDO®

Model	LEONARDO® DUAL 45
REF	SL980 + 1470nm45W
Wavelength	980 nm and 1470 nm
Power	max. 45 Watt (1470 nm/15 Watt + 980 nm/30 Watt) separately adjustable
Fiber diameter	≥ 360 µm
Aiming beam	532 nm and 635 nm, green 1 mW, red 4 mW, user controlled intensity
Treatment mode	CW, Pulse Mode, ELVeS® Signal, ELVeS® Segment, Derma Mode
Pulse duration /-break	0.01 – 60 sec / 0.01 – 60 sec
Power supply	110 – 240 VAC, 50 / 60 Hz, 450 VA
Batteries	–
Dimensions (H x W x D)	approx. 28 cm x 37 cm x 9 cm
Weight	approx. 8,5 kg

All laser sets incl. 3 safety goggles, foot switch, interlock connector, power cord and manual in a carrying case.

Litt fibers

REF	Product	PU [Packaging unit]	length [m]	ø fiber [mm]	Tip shape ablation
503200260	Litt fiber	5	10	1.85	 round
503200270	Litt2B fiber	5	10	1.85	 olive
503200280	Litt3B fiber	5	10	1.85	 olive

